



THERE'S ONE SMALL CAP THAT
SECURELY FITS BOTTLES AND VACUUM TUBES.

SampLok[®] Adapter Cap 3

Small, Safe, Simple Blood Culture Sample Collection

COMPACT SIZE generates reduced biohazard sharps waste

SECURE FIT makes it easy to handle while inserting,
filling and removing culture bottles and vacuum tubes

FLEXIBLE DESIGN allows sampling with a variety
of culture bottles and vacuum tubes



SampLok® Adapter Cap 3

Small, Safe, Simple Blood Culture Sample Collection

SampLok® Adapter Cap 3 is a versatile compact device for safe blood sample collection into culture bottles and vacuum tubes.

Fits securely over the blood culture bottle top.

Fits a variety of culture bottles and vacuum tubes.



PRODUCT FEATURES

The cap ensures the sampling needle is well positioned on the culture bottle septum and vacuum tube cap.

Small size reduces biohazard sharps waste.



Available bulk packaged or in shelf boxes.



Insert available for use with sample tubes and long neck culture bottles (optional).



Product Information

-  Not Made with Natural Rubber Latex
-  Single Use Only



Product Codes

A100760 Non-sterile, 600/case, 100/box, 6 boxes/case

A100761 Insert, non-sterile, 900/case, 150/box, 6 boxes/case

Can be attached to blood collection and winged needle sets.



HUMAN INSPIRED. SCIENTIFICALLY ENGINEERED.

ITL BioMedical

www.itlbiomedical.com

Sales Office: 1.888.411.2851

sales@itlbiomedical.com