



Vacuum Curettes


Fingertip precision

ITL BioMedical's range of Vacuum Curettes offer a choice of Straight, Curved or Flexible in multiple diameters.

Designed for patient comfort and clinician precision, the Vacuum Curettes are made of high quality materials to provide strength but minimise cervical dilation with thin wall design.

The textured band on the Curette allows for a visual and tactile indicator to gauge the depth of insertion.

**Partners in Healthcare Associated Blood
Stream Infection Prevention Products and
Strategies**



Patient Comfort and Clinical Precision

Vacuum Curettes

- Ergonomically designed intrauterine suction set with rotating hand piece increases functionality and precision with the inbuilt on/off slide control.
- Boasting a flute-shaped nozzle and a rounded tip, it seamlessly connects to the intrauterine suction tube.
- Extensive range in straight, curved and flexible available.
- Sterile and Single Use.
- Latex Free.
- Collection Set contains or presence of Bis(2-ethylhexyl) Phthalate (DEHP).
- All Curettes can be ordered in Boxes of 20ea.
- All Suction Sets can be ordered in Boxes of 10ea.



Product Range

Product Code	Diameter	Packaging
--------------	----------	-----------

Straight Rigid Vacuum Curettes

20/CRD07 S01	07mm	1/each
20/CRD08 S01	08mm	1/each
20/CRD09 S01	09mm	1/each
20/CRD10 S01	10mm	1/each
20/CRD11 S01	11mm	1/each
20/CRD12 S01	12mm	1/each

Also available in sizes 14 and 16mm.

Curved Rigid Vacuum Curettes

20/CRC07 S01	07mm	1/each
20/CRC08 S01	08mm	1/each
20/CRC09 S01	09mm	1/each
20/CRC10 S01	10mm	1/each
20/CRC11 S01	11mm	1/each
20/CRC12 S01	12mm	1/each

Flexible Vacuum Curettes

The Flexible ("Karman" style) cannulas have 2 openings allowing the cannula to collapse if pushed against the uterine fundus, minimising the risk of perforation.

20/CAN24 S01	04mm	1/each
20/CAN25 S01	05mm	1/each
20/CAN06 S01	06mm	1/each
20/CAN07 S01	07mm	1/each
20/CAN08 S01	08mm	1/each
20/CAN09 S01	09mm	1/each
20/CAN10 S01	10mm	1/each
20/CAN11 S01	11mm	1/each
20/CAN12 S01	12mm	1/each

Intrauterine Suction Set

The swivel handle contains a slip ring for fingertip vacuum control

20/SIVG1 S01	Collection Set, with handle, 2m length	1/each
20/SIVG2 S01	Collection Set, with handle, 3m length	1/each
20/SIVG3 S01	Collection Set, without handle, 2m length	1/each

For further information or samples, visit itlbiomedical.com or contact your local sales office.

Product Information



Made in Italy

www.itlbiomedical.com/symbolglossary

